

GIBCOBRLLIFE  TECHNOLOGIES™**SUPERSCRIPT™ II
RNase H⁻ Reverse Transcriptase**

Page 1 of 4

Cat. No. 18064-014**Lot No. _____ 10,000 units; 200 U/μl****Exp. Date: ____ . Store at -20°C (not frost-free).**

SUPERSCRIPT II RNase H⁻ Reverse Transcriptase (U.S. Patent 5,244,797) is purified to near homogeneity from *E. coli* containing the *pol* gene of Moloney Murine Leukemia Virus (1,2). The enzyme is used to synthesize first strand cDNA and will generally give higher yields of cDNA and more full length product than other reverse transcriptases.

Components:

18064-014	SUPERSCRIPT II	Lot No.
Y00146	5X First Strand Buffer	Lot No.
Y00147	0.1 M DTT	Lot No.

Unit Definition:

One unit incorporates 1 nmole of dTTP into acid-precipitable material in 10 min at 37°C using poly(A)•oligo(dT)₂₅ as template-primer (3).

Storage Buffer:

20 mM Tris-HCl (pH 7.5)
100 mM NaCl
0.1 mM EDTA
1 mM DTT
0.01% (v/v) NP-40
50% (v/v) glycerol

Additional Components Provided:

The 5X First Strand Buffer [250 mM Tris-HCl (pH 8.3 at room temperature), 375 mM KCl, 15 mM MgCl₂] and 0.1 M DTT are provided. Store both solutions at -20°C. Thaw the solutions at room temperature just prior to use and refreeze immediately.

Doc. Rev.: 101397

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

For technical questions about this product, call the Life Technologies TECH-LINE™ (800) 828-6686.

Quality Control Assays:

This product has passed the following quality control assays:

SDS-polyacrylamide gel analysis for purity; functional absence of endodeoxyribonuclease, 3' and 5' exodeoxyribonuclease, and ribonuclease activities; yield and length of cDNA product.

Note: Limited Label License

This product is the subject of claims in US Patent 5,405,776 and 5,244,797, owned by Life Technologies, Inc. Purchase of this product conveys to the buyer only the non-transferable right under these patents to use the product in research conducted by the buyer. The buyer cannot sell or otherwise transfer this product to a third party. Life Technologies, Inc. reserves all other rights, and this product may not be used in any manner other than provided herein. In particular, without prior written agreement between the buyer and Life Technologies, Inc., no rights under these above patents are conveyed to the buyer to use the product in manufacturing, or to resell the product or components of the product as a standalone reagent or in a kit, whether or not such standalone reagent is resold or such kit is resold for use in research. If the purchaser is not willing to accept the limitations of this Label License, Life Technologies, Inc. is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact the Director of Corporate Development, 9800 Medical Center Drive, Rockville, MD 20849. Phone (301) 610-8000. Fax (301) 610-8383.

First Strand cDNA Synthesis Using SUPERSCRIPT II for RT-PCR:

A 20- μ l reaction volume can be used for 1-5 μ g of total RNA or 50-500 ng of mRNA. Add the following components to a nuclease-free microcentrifuge tube:

1 μ l Oligo (dT)₁₂₋₁₈ (500 μ g/ml)^a

1-5 μ g total RNA

Sterile, distilled water to 12 μ l

Heat mixture to 70°C for 10 min and quick chill on ice. Collect the contents of the tube by brief centrifugation and add:

4 μ l 5X First Strand Buffer

2 μ l 0.1 M DTT

1 μ l 10 mM dNTP Mix (10 mM each dATP, dGTP, dCTP and dTTP at neutral pH)

Mix contents of the tube gently and incubate at 42°C for 2 min. Add 1 μ l (200 units) of SUPERSCRIPT II, mix by pipetting gently up and down. Incubate 50 min at 42°C. Inactivate the reaction by heating at 70°C for 15 min. The cDNA can now be used as a template for amplification in PCR. However, amplification of some PCR targets (those > 1 kb) may require the removal of RNA complementary to the cDNA. To remove RNA complementary to the cDNA, add 1 μ l (2 units) of *E. coli* RNase H and incubate 37°C for 20 min.

^a Alternatively 50-250 ng of random primers or 2 pmole of a gene specific primer may be used. Use of random primers requires incubation at 25°C for 10 min before the 42°C incubation.

PCR Reaction

Use only 10% of the first strand reaction for PCR. Adding larger amounts of the first strand reaction may not increase amplification and may result in decreased amounts of PCR product.

1. Add the following to a PCR reaction tube for a final reaction volume of 100 μ l:

10 μ l	10X PCR Buffer [200 mM Tris-HCl (pH 8.4), 500 mM KCl]
3 μ l	50 mM MgCl ₂ ^b
2 μ l	10 mM dNTP Mix
1 μ l	amplification primer 1 (10 μ M)
1 μ l	amplification primer 2 (10 μ M)
1 μ l	Taq DNA polymerase (5 U/ μ l)
2 μ l	cDNA (from first strand reaction)
80 μ l	autoclaved, distilled water

2. Mix gently and layer 2 drops (~100 μ l) of silicone oil over the reaction. (*Note: the addition of silicone oil is unnecessary in thermal cyclers equipped with a heated lid.*)
3. Heat reaction to 94°C for 3 min to denature.
4. Perform 15 to 40 cycles of PCR. Annealing and extension conditions are primer and template dependent and must be determined empirically.

^b Optimal concentration of MgCl₂ needs to be determined empirically for each template primer pair.

References:

1. Kotewicz, M.L., D'Alessio, J.M., Driftmier, K.M., Blodgett, K.P., and Gerard, G.F. (1985) *Gene* 35, 249.
2. Gerard, G.F., D'Alessio, J.M., Kotewicz, M.L., and Noon, M.C. (1986) *DNA* 5, 271.
3. Houts, G.E., Miyagi, M., Ellis, C., Beard, A., and Beard, J.W. (1979) *J. Virol.* 29, 517.